

PATROLS

Advanced Tools for NanoSafety Testing

This project has received funding
from the European Union's Horizon
2020 research and innovation
programme under grant agreement
No 760813



PATROLS Standard Operating Procedures (SOP)

Guidance Document for DQ₁₂ and TiO₂ Aerosolization using VITROCELL[®] Cloud System

SOP_PATROLS_Cloud_Aerosolization

**This is a SOP recommended for
external use by PATROLS**

Adapted from the NanoImpactNet SOP, Clift *et al* (Deliverable 5.4 under the European Commission's 7th Framework Programme, Grant Agreement 218539).

This is an Open Access document distributed under the terms of the Creative Commons Attribution-NonCommercial-ShareAlike 4.0 International License. To view a copy of this license, visit <http://creativecommons.org/licenses/by-nc-sa/4.0/>, which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited.

Authored by:

Hana Barosova, Anne Bannuscher, Barbara Drasler, Barbara Rothen-Rutishauser

Adolphe Merkle Institute (AMI), University of Fribourg, Switzerland

Reviewed by:

PATROLS WP3 and external members

Document History:

Version	Approval Date	Description of the change	Author(s) of change
1.0	26.07.2019	Initial Document	Hana Barosova
1.1	04.03.2020	Additional detail added by AMI	Anne Bannuscher
1.2	20.5.2020	Additional detail added by AMI	Barbara Drasler
1.3	25.06.2020	Additional details added by AMI	Anne Bannuscher, Barbara Drasler
1.4	29.06.2020	Minor corrections by RIVM	Hedwig Braakhuis
1.4	01.07.2020	Final coments implemented	Barbara Drasler
1.4	08.07.2020	Version distributed to WP3 members	Anne Bannuscher
1.5	14.07.2020	Additional details added	Otmar Schmid, Yaobo Ding, Anne Bannuscher, Barbara Drasler
2.0	21.07.2020	Corrected version distributed to the WP3 partners	Anne Bannuscher
2.1	21.09.2020	Corrected version distributed to the WP3 partners	Anne Bannuscher, Barbara Drasler
2.2		All commets from WP3 members intergrated and uploaded to server	
2.3	15.10.2022	Adaptation to a general-use SOP: Clarifications and edits were added	Otmar Schmid, Barbara Rothen-Rutishauser

Table of Contents

Contents

1	INTRODUCTION	4
1.1	SCOPE AND LIMITS OF THE PROTOCOL	4
2	TERMS AND DEFINITIONS	5
3	ABBREVIATIONS	6
4	PRINCIPLE OF THE METHOD	6
5	DESCRIPTION OF THE METHOD	7
5.1	TEST SYSTEM USED	7
5.2	CHEMICALS, NANOMATERIALS AND REAGENTS USED	7
5.3	APPARATUS AND EQUIPMENT USED	7
5.4	REPORTING OF PROTECTED ELEMENTS.....	8
5.5	HEALTH AND SAFETY PRECAUTIONS.....	8
5.6	PARTICLE SUSPENSION PREPARATION	9
6	PROCEDURE	10
6.1	PRIOR TO AEROSOLISATION	10
6.2	AEROSOLISATION: GENERAL PROCEDURE	11
6.3	SPECIFIC REMARKS FOR PARTICLE AEROSOLIZATION.....	13
6.4	CLEANING.....	16
7	DQ₁₂ CONCENTRATIONS APPLIED AND DEPOSITION MEASURED	17
8	TIO₂ NM-105 CONCENTRATIONS APPLIED AND DEPOSITION MEASURED	18
9	TEM IMAGES	19
10	FLUORESCHEIN ASSAY REMARKS	19
11	QUALITY CONTROL & ACCEPTANCE CRITERIA	20
12	DATA ANALYSIS AND REPORTING OF DATA	20
13	PUBLICATIONS	20
14	REFERENCES	21

1 Introduction

Due to the constant increase in their production, exposure to engineered nanomaterials (ENM) poses an inevitable health risk to both humans and the environment through long-term, repetitive, low-dose exposures. The majority of the literature however, focuses on short-term, high-dose exposures. Hazard assessment of ENM, when applying alternative testing strategies to *in vivo* research, has previously engaged exposures under submerged conditions. Such exposures poorly mimic human exposures, and possess several limitations, such as issue with measuring a delivered dose, and interaction with exposure medium. Therefore, air-liquid interface (ALI) exposure chambers allowing for single droplet deposition of ENM aerosols onto cell surface were developed. VITROCELL[®] Cloud system equipped with vibrating mesh nebulizer (Aeroneb Pro/Lab, Aerogen Inc., Ireland), and Quartz crystal microbalance (QCM) for measuring online deposition of non-/poorly-soluble (in water) particles such as titanium dioxide (TiO₂) NM-105, crystalline SiO₂ (DQ₁₂ quartz) and substances such as sodium fluorescein salt, is referred in this SOP. A QCM coupled to Cloud12 present a highly sensitive, accurate device for (quasi-) real-time dosimetry of the cell-delivered particle dose in ALI cell exposure experiments. In addition, for a verification on the QCM's accuracy a fluorescein-based method is included as quality control to prevent measurement errors due to damaged QCM crystals or electronic failurs. Detailed descriptions of the test system can be found in references [1-3].

1.1 Scope and limits of the protocol

This SOP was established with the intention to be used for developing an exposure protocol for a lung cell culture model using a VITROCELL[®] Cloud exposure system which can provide a physiologically relevant assesement of the hazards associated with ENM exposures over both an acute and chronic, repeated dose regime. This SOP provides instructions on how to aerosolize non- or poorly-soluble (in water) particles like DQ₁₂ or ENM such as TiO₂ NM-105 and water-soluble/-dissolved (molecular) substances such as fluorescein sodium salt during a interlaboratory test within the H2020 funded PATROLS project (Physiologically Anchored Tools for Realistic nanOmateriAL hazard aSsessment, No. 760813); for handling the VITROCELL[®] Cloud

system in general please refer to the instructions provided by VITROCELL® SYSTEMS (Germany).

Limitations:

The VITROCELL® Cloud system requires preparation of a stable ENM (or particle) suspension. For this, the NanoReg protocol was used: “Protocol for producing reproducible suspensions of manufactured nanomaterials in environmental exposure media” (https://www.anses.fr/en/system/files/nanogenotox_deliverable_5.pdf). Details on the suspension preparation are given in section 5.6. Depending on the properties of the ENM or particles, the NanoReg protocol may not allow generating a stable ENM/particle suspension (i.e., not all of the ENM/particles can be dispersed uniformly in aqueous liquid without ENM/particle settling for a long enough period to transfer the suspension into aerosol with the VITROCELL® Cloud system). Then, the NanoReg protocol has to be adapted or a different protocol has to be used.

In general, aerosolization of ENMs/particles with the VITROCELL® Cloud system is limited by the pore size of its vibrating mesh nebulizer. For the typically used Aeroneb Pro/Lab nebulizers this typically requires that the ENM/particles and agglomerates thereof should be smaller than the diameter of the pores of the vibrating mesh nebulizer (here about 5 µm) in at least one dimension. Thus, nebulization of, e.g., fiber-like particles larger than 5 µm is possible, if their cross-sectional diameter is smaller than ca. 5 µm.

2 Terms and Definitions

Nanomaterial

Material with at least one external dimension in the *nanoscale* or having internal structure or surface structure in the nanoscale.

Note 1 to entry: This generic term is inclusive of *nano-object* and *nanostructured material*.

Note 2 to entry: In the context of this SOP, the material is non- or poorly-soluble in water, i.e., it maintains its integrity when dispersed in aqueous liquid.

[SOURCE: ISO/TS 80004-1: 2016, definition 2.4]

Engineered nanomaterial

Nanomaterial designed for specific purpose or function

[SOURCE: ISO/TS 80004-1: 2016, definition 2.8]

Particle

Minute piece of matter with defined physical boundaries, which is not or only poorly soluble in aqueous liquids.

Note 1 to entry: A physical boundary can also be described as an interface.

Note 2 to entry: A particle can move as a unit.

Note 3 to entry: This general particle definition applies to *nano-objects*.

[SOURCE: ISO 26824:2013, 1.1]

Substance

Single chemical element or compound, or a complex structure of compounds.

3 Abbreviations

ALI – Air-liquid interface

BSA – Bovine Serum Albumine

DQ12 - Dörnruper quartz 12

ENM – Engineered Nanomaterials

QCM – Quartz crystal microbalance

TEM – Transmission electron microscopy

4 Principle of the Method

This method aims to standardise the aerosolization of ENM using VITROCELL® Cloud system.

This protocol will be broken into key stages:

1. Particle suspension preparation
2. Preparing the aerosolization device
3. Aerosolization
4. Cleaning

5 Description of the Method

5.1 Test system used

This SOP should be carried out under laboratory based conditions, with all work following safe handling of ENM.

- VITROCELL[®] Cloud12 system (VITROCELL[®] SYSTEMS, Germany) equipped with
 - Quartz crystal microbalance (QCM)
 - Aeroneb[®] Lab nebulizer (4-6 µm pore size, Aerogen, Ireland)
 - This protocol refers to Cloud 12, which is designed for 12-well inserts.
- For further information:
<https://www.vitrocell.com/inhalation-toxicology/exposure-systems/vitrocell-cloud-system/vitrocell-cloud-12>

5.2 Chemicals, Nanomaterials and Reagents used

- Ultrapure water
- 70% ethanol or a similar cleaning agent (not Aceton or Sodium-Hydroxide; usage of those will cause unrepairable damage to the system or parts of it)
- Isotonic (0.9%) NaCl solution (90 mg NaCl/ mL water)
- Silica particles - DQ₁₂ (Institute of Occupational Medicine, United Kingdom)
- TiO₂ NM-105 (Fraunhofer Institute for Molecular Biology and Applied Ecology (IME), Germany)
- Fluorescein sodium salt (CAS 518-47-8)
- Dulbecco's Phosphate Buffered Saline (DPBS Gibco, 14190144)

5.3 Apparatus and equipment used

- 10 µL, 200 µL and 1000 µL pipette tips

- 50 mL centrifuge tubes
- Glas vial e.g. 20 mL Scint-Burk glass pp-lock+Alu-foil vials (WHEA986581; Wheaton Industries Inc.)
- Sonicator (Branson 550)
- VITROCELL® Cloud12
- Aeroneb Pro/Lab vibrating mesh nebulizers (4-6 µm pore size) (Aerogen Inc, Galway, Ireland)
- Cell culture inserts (e.g. Corning Transwell inserts for 12-wells plates, 0.4 µm pore size (Cat.no 3460))
- TEM grids (EMS, CF300-Cu)
- TEM grid holder for 12-wells (VITROCELL® SYSTEMS, VC5036)
- Water bath sonifier (e.g. Elmasonic P30H cleaning unit, 100 W, 37 kHz, 30% (Elma Schmidbauer GmbH, Germany))
- Microbalance (integrated in VITROCELL® Cloud)
- Wipers e.g. Kimberly-Clark Professional™ Kimtech Science™ Precision Wipes™ Tissue Wipers (19-063-099) or similar lint-free cloth

5.4 Reporting of protected elements

To the best of our knowledge, this SOP does not have any associated patent restrictions, specific licenses, material transfer agreements or commercial purchase requirements required to perform the protocol described.

5.5 Health and safety precautions

Standard health and safety precautions associated with working within a laboratory environment and handling ENM, as described by the European Agency for Safety and Health at Work (<https://osha.europa.eu/en/safety-and-health-legislation/european-guidelines>), should be adopted when conducting this SOP. In addition, all health and

safety precautions outlined in the MSDS data sheets associated with the specific chemicals required must also be followed. DQ₁₂ is a non-nanosized particle which can induce inflammation and fibrosis in humans upon inhalation. Therefore, suitable protective equipment such as a lab coat, goggles, gloves and a mask must be worn.

5.6 Particle suspension preparation

Weighing of DQ₁₂ and ENM such as TiO₂ NM-105 should be performed in a ventilated weighing box, glove box or fume hood designed for sensitive weighing. Ensure appropriate safety clothing, including the use of double gloves. DQ₁₂ and ENM suspensions have to be prepared in ultrapure water (**without** BSA, which is different to the Nanogenotox protocol) with EtOH pre-wetting. Suspensions have to be prepared in a volume range between 4 mL and 10 mL, following the Nanogenotox protocol (https://www.anses.fr/en/system/files/nanogenotox_deliverable_5.pdf).

1. Using a microbalance, weight the respective amount of DQ₁₂/ENM in a glass vial corresponding to a minimum of 4 mL to a maximum of 10 mL (e.g. 10.24 mg corresponding to the final volume of 4 mL).
2. Pre-wet the particles with ethanol (e.g. 20 µL corresponding to 10.24 mg and 4 mL or 30 µL for 15.36 mg and 6 mL), add the corresponding amount of ultrapure water (in 4 steps) and mix it carefully by vortexing to obtain a 2.56 mg/mL stock solution.
3. Sonicate the suspension based on the Nanoreg protocol, please calibrate your sonifier respectively as described in “NANoREG SOP for probe-sonicator calibration of delivered acoustic power and de-agglomeration efficiency for in vitro and in vivo toxicological testing” (4).
4. Carefully remove the sonicator tip, close the vial and clean the device.
5. Store the aliquots at 4°C until use (suspensions need to be prepared freshly on the day of experiment).
6. Prepare a 1 mg/mL dilution from the stock solution (e.g. 390.6 µL of the 2.56 mg/mL stock and 609.4 µL of ultrapure water). Dilute this stock further to the desired concentrations (250 µg/mL or 500 µg/mL).

6 Procedure

6.1 Prior to aerosolisation

1. Perform all the steps in a laminar flow hood.

Note: Figure 1 shows the VITROCELL[®] Cloud system with all of its parts.

2. Assemble the QCM, integrate it into one of the wells of the VITROCELL[®] Cloud system and connect it with the oscillator.

Note: A YouTube video is available

<https://www.youtube.com/watch?v=HE0AGAvabck>, accessed 23 July 2022.

3. Turn on the VITROCELL[®] Cloud heating unit and wait until the temperature reaches 37°C (typically ca. 30-45 minutes).
4. Warm the DQ₁₂ or ENM suspension to room temperature and shortly sonicate (ultrasound bath) for 10 min to redisperse the sedimented/agglomerated particles or use them directly after the preparation in the section 5.6. For the purpose of this SOP, each particle suspension has to be prepared freshly before usage.




Figure 1: The VITROCELL® Cloud 12 system.

6.2 Aerosolisation: general procedure

The general procedure for performing the aerosolisation is described below (see the schematic depiction in Figure 2).

1. Rinse the reservoir of nebulizer with ultrapure water and dry it with lint-free cloths. Place the nebulizer on top of the exposure chamber.
2. Prepare a solution of ultrapure water or DQ₁₂/ENM suspension containing 1% of isotonic NaCl solution (v/v), e.g. 2.5 μL in 250 μL and pipette 200 μL to the nebulizer (make sure the liquid covers the nebulizer's mesh).
3. Start the data acquisition on the QCM and wait until the signal is stable (typically ca. 30 s, but can be longer), set the signal to zero and wait another 30 s to verify stability of the zero point signal (the signal should be in the range of ± 25 ng/cm^2). If this is not the case, set the QCM signal to zero again and wait for another 30 s to verify stability of the signal.

Note: It is recommended to zero the signal once again after starting the data acquisition.

4. Measure the signal for 1 minute before starting the nebulization (next step).
5. Start the nebulization process by plugging in the nebulizer in the driving unit and at the same time start the timer ($t = 0$ min). Carefully measure the time needed for nebulization by watching either the cloud emitted from the nebulizer or the decrease/vanishing of particle suspension in the reservoir of the nebulizer.

Note: The acceptable time range for nebulizing 200 μL of water/suspension is between 15-60 s (0.2 – 0.8 mL/min).

Note: The output rate should not exceed 10% deviation from the original output rate of the individual nebulizer and this should be measured prior to all nebulizations. Prolonged nebulization time can be a sign of a blocked nebulizer requiring thorough cleaning (ultrabath sonication, see section 6.3) or replacement.

6. Wait for 6 min to ensure a complete settlement of the test substance onto the inserts and QCM. A typical time course of the QCM signal is seen in Figure 2.

Note: If the QCM signal is not dropping within this time frame, extend this waiting period to 10 min and follow this procedure for all nebulizations

7. Lift the exposure top for 1 min allowing the QCM to dry.

Note: Complete drying of the QCM is typically associated with a step-change in QCM signal (either increase or decrease (as in Figure 2) of signal depending on the deposited particle dose). In case of conduction of further nebulizations, the inner walls of the exposure chamber need to be wiped with 70% ethanol.

8. After this 1 min place the exposure top back onto the device.

Note: The reason behind covering the device is to minimize the effects of ventilation and temperature fluctuation on the measured QCM values.

9. Wait for 3 min and stop the data acquisition.

10. Lift the exposure top and clean the QCM: use a Q-tip wrapped in lint-free wipers for cleaning thoroughly the crystal.

Note: A YouTube video is available

<https://www.youtube.com/watch?v=HE0AGAvabck>, accessed 23 July 2022

Note: Vitrocell recommends to change crystal after each nebulization

11. Data analysis should be performed on the mean of the last 30 seconds of the QCM signal (here: last 30 data points).

Note: The accuracy of Cloud12 QCM is $3.8\% \pm 1.0\%$ [3]




Figure 2: Typical time course of QCM signal and the time periods applied.

6.3 Specific remarks for particle aerosolization

NOTE: Please see sections 7 and 8 for recommended concentrations of particle suspensions and the respective deposited values

1. Prior to each aerosol exposure run, 200 μL ultrapure water containing 1% isotonic NaCl solution should be nebulized (into the hood – not into the VITROCELL[®] Cloud system) for removal of residual ENM from the nebulizer. Dry the nebulizer using the lint-free wipers (reservoir and from the bottom side – gently touching the vibrating mesh) in order to remove residual liquid.
2. At the beginning of a working day, measure the deposited dose of 200 μL of ultrapure water with 1% (v/v) isotonic NaCl solution (blank value).

Note: For 200 μL of water with 1% isotonic NaCl solution a deposited dose of around 100 ng/cm^2 is expected. However, due to the proximity of this value to the detection limit (169 ng/cm^2 [3]), the values could also be negative or significantly larger than 100 ng/cm^2 (up to 250 ng/cm^2). If the values are above the threshold of 250 ng/cm^2 , this could be a reflection either of a dirty nebulizer, i.e. there are remnants of material from previous nebulizations, or of trouble with the QCM signal. Please note, that the blank value can be strongly biased by small shifts in the zero point signal of the QCM. Thus, it is recommended to set the zero point of the QCM again to zero and measure the deposited dose of 200

μL of ultrapure water with 1% (v/v) isotonic NaCl solution (blank value) again. If the QCM signal is still above the threshold of 250 ng/cm^2 , carefully clean the nebulizer (see section 6.3) and redo nebulization of $200 \mu\text{L}$ of water with 1% isotonic NaCl solution. If the blank value of the QCM still remains above 250 ng/cm^2 take the QCM apart and clean all its parts (see the next step and the video: <https://www.youtube.com/watch?v=HE0AGAvabck>, accessed 23 July 2022).

3. QCM cleaning between subsequent nebulizations: Wipe the entire QCM crystal carefully with water on a Q tip, also wipe the edges and every surface that came in contact with nebulized material very carefully (but not too harshly). Repeat the procedure with 70% ethanol. Then wait at least 5 min to ensure complete drying of the crystal as reflected by a stable QCM signal.

Note: The quartz crystal of the QCM needs to be replaced when scratches are observed on the surface and/or if a stable signal cannot be achieved after the crystal and base module have been heated up to 37°C with the exposure top placed on the base module.

4. Place one TEM grid in the TEM grid holder (provided by VITROCELL[®]) for each technical repetition, and place the TEM grid holder in one well of the base module.
5. Prepare $\text{DQ}_{12}/\text{ENM}$ suspension in water with the desired concentration (see section 5.6).
6. Mix $250 \mu\text{L}$ of $\text{DQ}_{12}/\text{ENM}$ suspension with $2.5 \mu\text{L}$ of isotonic NaCl solution in an Eppendorf tube; for the concentrations see sections 7 and 8.
7. Pipette $200 \mu\text{L}$ of the suspension into the nebulizer reservoir. Keep the nebulizer upright while placing it into the circular gap in the exposure top. Make sure, that the nebulizer is placed straight, and sits tightly. Then connect it and start nebulization.

Note: Follow the data acquisition, nebulization and timing procedure as indicated in the steps 5-10 of the general procedure and in the Figure 2.

8. Transfer the TEM grid onto a dust-protected surface allocated for drying, e.g. a parafilm sheet, and let it dry at least over night before transferring it into a TEM grid box.
9. Wash the nebulizer with 200 μ L ultrapure water containing 1% (v/v) isotonic NaCl solution. Measure output rate of the nebulizer after each nebulization (should not exceed 10% deviation from the original output rate).

Note: If the time is exceeding, put the nebulizer in a beaker with ultrapure water or 70% ethanol and sonicate it in an ultrasonic bath for 5-10 min. Then, measure the output rate again. If nebulizing time is not improving, clean more thoroughly (see section 6.3) or replace the nebulizer.

If no cell exposure work is involved, cleaning in water is sufficient. Otherwise, an additional step of ultrabath sonication in 70% ethanol (e.g. 10 min) for sterilization is required.

10. Wipe the nebulizer reservoir and the bottom part of the nebulizer in order to remove remaining liquid (use lint-free wipes).
11. Continue with the step 6, i.e. nebulize the same particle suspensions for technical repetitions or other concentrations.

Note: Clean the nebulizer carefully if different substances i.e. particles are to be nebulized, see section 6.4 Cleaning between successive runs of the same substance can prevent blocking but is optional.

12. For each particle 3 technical replicates have to be performed for the selected concentration.
13. At the end of a working day and/or one type of material, cleaning steps need to be performed, see section 6.4 Cleaning.
14. For subsequent nebulizations of the same particle solution concentration: in order to achieve a higher deposition subsequent nebulizations are possible, but to a maximum of 4. After the first nebulization, check the output rate of the nebulizer. If this is within 10% of the initial output rate, nebulize 1x 200 μ l ultrapure water with salt (not into the cloud) and continue with the second

nebulization. If not, put the nebulizer in a beaker with ultrapure water or 70% ethanol and sonicate it in an ultrasonic bath for 5-10 min and proceed. **Do not stop the data acquisition and do not clean the QCM in between** (possible to a maximum deposition of $\sim 1.5 \mu\text{g}/\text{cm}^2$). The measured deposited dose should be approximately twice the first measured dose. If not, the nebulizer may be blocked and must be cleaned carefully, see section 6.3.

6.4 *Cleaning*

15. Nebulize at least 3x 200 μL of ultrapure water containing 1% (v/v) isotonic NaCl solution.
16. Put the nebulizer in a beaker with ultrapure water in water bath sonicator (e.g. Elmasonic P30H cleaning unit, 100 W, 37 kHz, 30% (Elma Schmidbauer GmbH, Germany)) for 10 min.
17. Rinse the nebulizer with ultrapure water.
18. Put the nebulizer in a beaker with 70% ethanol and sonicate in water bath sonicator (e.g. for 10 min).
19. Clean VITROCELL[®] base module and cover top with 70% ethanol.
20. QCM Crystal: Used, dirty crystals can be collected, cleaned and re-used, if no scratches on the surfaces are visible. For cleaning put them in ethanol for 0.5 min to 5 min (depends on the degree of impurity/dirt) and wipe them carefully afterwards. Check their performance afterwards and trash them if signal oscillates too strongly, loss of signal is present or scratches on the surface are visible.
21. VITROCELL[®] base module and exposure top can be washed in dishwasher (recommended). Do NOT clean with acetone.
22. VITROCELL[®] base module can be autoclaved (recommended on regular basis, *i.e.*, every week, but not necessary after each experiment) following intervals and temperatures listed, *i.e.*, minimum warm-up of 30 min, autoclaving at 121°C for a minimum of 20 min, minimum cool-down time in the autoclave of 30 min.

7 DQ₁₂ concentrations applied and deposition measured

For the SOP optimization purposes, 3 individual nebulizations of the concentration 500 µg/mL of DQ₁₂ are required. Please provide also the deposition value of ultrapure water containing 1 % (v/v) isotonic NaCl solution (this needs to be done just before nebulizing DQ₁₂ using the same nebulizer and crystal).

Table 1: Deposited values of nebulized DQ₁₂ in a concentration range. N=6

Concentration (µg/mL)	Deposition (ng/cm ²) ± SD
water	44.9 ± 51
125	185.2 ± 19
250	258.3 ± 19
500	516.0 ± 32

8 TiO₂ NM-105 concentrations applied and deposition measured

For the SOP optimization purposes, 3 individual nebulizations of the concentration 500 µg/mL of TiO₂ NM-105 are required. Please provide also the deposition value of ultrapure water containing 1 % (v/v) isotonic NaCl (this needs to be done just before nebulizing TiO₂ NM-105 by using the same nebulizer and crystal).


Note: If the concentration 500 µg/mL of TiO₂ NM-105 is causing blocking of the nebulizer, the nebulizer needs to be thoroughly cleaned and the concentration 250 µg/mL shall be measured.

Table 2: Deposited values of nebulized TiO₂ NM-105 in a concentration. N=6


Concentration (µg/mL)	Deposition (ng/cm ²) ± SD
water	50.8 ± 36
125	192.5 ± 46
250	217.6 ± 42
500	394.8 ± 121

9 TEM images


DQ₁₂



125




250




500 µg/ml


TiO₂ NM-105



125 µg/ml



250 µg/ml



500 µg/ml

10 Fluorescein assay remarks

The fluorescein assay was performed by each partner following the guidance document from VITROCELL® Systems (The document is available for VITROCELL® customers. Contact: support@vitrocell.com). The purpose of the assay is to assess the deposition factor of the system which enables to calculate the delivered dose of aerosols to the cells.

Important remarks to the SOP:

- Use 0.4 µm pores in the transwell insert membrane and AVOID any contact of inserts with otherwise the chemical can migrate into PBS and DF values are incorrect.
- Make sure that the device is heated to 37 °C.

- Perform the experiments in a laminar flow hood (the same place and position as for the DQ₁₂/ENM nebulization).
- In case the instrument's (fluorimeter's) detection limit allows for, use 15 µg/ml fluorescein sodium salt in DPBS as described in the SOP; otherwise use 30 µg/mL.
- 3 repetitions need to be performed

A brief description of the procedure (please refer to the guidance document from VITROCELL® Systems for the fluorescein assay for a detailed description):

- 1.) Prepare fluorescein sodium salt stock solution (i.e. 15 µg/mL).
- 2.) Heat the device to 37°C including QCM.
- 3.) Fill the wells of the basal module with ultrapure water (3 mL) and the top chamber of the inserts with DPBS (0.3 mL).
- 4.) Perform fluorescein sodium salt nebulization. Try to keep the time between filling liquids and aerosolisation as short as possible. Otherwise, signal losses can follow.
- 5.) Follow exactly the same data acquisition pattern as presented in section 6.2, i.e. 1 min QCM signal acquisition before the start of nebulization, wait for 6 min to ensure the aerosol is completely settled, 0.5 min of lifting the exposure top, followed by covering it back for additional 2.5 min to stabilize the signal.

11 Quality control & acceptance criteria

The deposited particles should be checked using TEM after trial exposures, before starting cell culture experiment, and then several times during the experiment.

12 Data Analysis and Reporting of Data

Not applicable for this current SOP.

13 Publications

Not applicable for this current SOP.

14 References

- [1] Lenz, A.G., Karg, E., Lentner, B. et al. A dose-controlled system for air-liquid interface cell exposure and application to zinc oxide nanoparticles. *Part Fibre Toxicol* 6, 32 (2009). <https://doi.org/10.1186/1743-8977-6-32>
- [2] Lenz AG, Stoeger T, Cei D, Schmidmeir M, Semren N, Burgstaller G, Lentner B, Eickelberg O, Meiners S, Schmid O. Efficient bioactive delivery of aerosolized drugs to human pulmonary epithelial cells cultured in air-liquid interface conditions. *Am J Respir Cell Mol Biol*. 2014 Oct;51(4):526-35. doi: 10.1165/rcmb.2013-0479OC.
- [3] Ding, Y., Weindl, P., Lenz, A. *et al.* Quartz crystal microbalances (QCM) are suitable for real-time dosimetry in nanotoxicological studies using VITROCELL®Cloud cell exposure systems. *Part Fibre Toxicol* 17, 44 (2020). <https://doi.org/10.1186/s12989-020-00376-w>
- [4] Booth A, Jensen KA. NANoREG D4.12 SOP Probe Sonicator Calibration for Ecotoxicological Testing. Available online: <https://www.rivm.nl/en/documenten/nanoreg-d412-sop-probe-sonicator-calibration-for-ecotoxicological-testing> (accessed on 18 July 2022).